

Contactless full scale semiconductor wafer chemical composition scanner

## EURO-X-SCAN

**CONTACTLESS FULL SCALE SEMICONDUCTOR WAFER CHEMICAL COMPOSITION IMAGING BY CAPILLARY FOCUSED X-RAY MICRO BEAM FLUORESCENCE.**

X-ray fluorescence is a widely used technique for chemical composition analysis.

The fluorescence is generally excited through a shielded X-ray beam or an electronic microscope, however the shielding reduces drastically the beam intensity and the electronic microscope needs vacuum operation. During the growth of compound semiconductor crystals, both deviations from composition are observed, resulting from microscopic changes (for example due to the segregation coefficients) as well as microscopic scale (non-uniformity, precipitates, etc...).

These composition changes have drastic effects on the characteristics of the materials and consequently on the performance of the devices.

# EURO-X-SCAN

## NO DIRECT CONTACT WITH THE WAFERS


**THE STANDARD VERSION OF THE INSTRUMENT** is able to measure with a special resolution of typically 20µm the X-ray excited fluorescence lines of the chemicals, without any physical contact on a full scale wafers by scanning.

But, besides semiconductor wafers, the **EURO-X SCAN** can be used to **analyse** with a high local resolution **any solid surface**.


**THE PRINCIPLE OF THE MEASUREMENT** is to excite locally the characteristic X-ray lines by a beam of X-rays given by a tube and having been focused to diameters in the range of 20µm through a special capillary.

## THE INSTRUMENT CONSISTS OF THE FOLLOWING PARTS:

- an X-ray generator coupled to
- a capillary tube, based on critical angle full reflectivity of the photons
- a computer controlled X-Y displacement system
- a high resolution dispersive semiconductor spectrometer, X-ray spectroscopy system operating at RT
- a PC based software




Principle of the method




Technical shematic of the system

# EURO-X-SCAN



## THE EURO-X-SCAN MEASUREMENT SYSTEM


## THE EURO-X-SCAN MEASUREMENT CHAMBER



# EURO-X-SCAN

## RESULTS

To illustrate the possibilities of the instrument, CZT (CdZnTe) wafers have been analysed, cut both perpendicular and parallel to the growth direction.



Control & measurement window

## 2D (Zn) MAPPING OF CdZnTe WAFERS to visualize changes in the Zn concentration over the whole ingot

